



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Transport and hospital trolley (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- **New mattress upholstery with medical certificate in green,**
- 2-segment transport trolley for free movement of patients around medical facilities and clinics,
- Foldable side rails,
- Built-in bumpers,
- **Manual adjustment:**
 - Gradual adjustment of the headrest backrest in the range of: 48°,
 - Free adjustment of the bed height in the range: 61 – 81 cm,
- **Mobile transport trolley on 4 wheels (2 braking),**
- **Included in the set:**
 - Patient mattress dimensions: 60 x 178 cm,
- Overall dimensions: 66 x 190 x 61/81 cm,
- Weight: 40 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for domestic market (Poland),**
 - **3 months to international market,**
 - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

