



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



LOGIQ P5 ultrasound machine (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Japanese production, (2014),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 750 VA,
GE Logiq P5
is a high-performance ultrasound machine with shared services. It supports a variety of applications, including OB/GYN and cardiology, It can perform 4D obstetric scans (optional), support surgical procedures and perform a variety of cardiac scans. GE Logiq P5's premium features include SRI (speckle reduction imaging), CE (coded excitation), ASO (automatic spectrum optimization), Crossbeam technology and more,
GE Logiq P5
is designed with a simple, ergonomic layout that enables fast workflow and productivity. It features a 15-inch TFT monitor with a fully movable arm for better viewing angles for sonographers and patients,
Specifications:
15-inch LCD monitor on an articulating arm,
Full-size keyboard with user-defined keys,
Automatic Optimization (AO),
Speckle Reduction Imaging (SRI),
3D/4D imaging (appropriate head required),
TrueScan architecture,
True Access raw data processing,
SmartScan,
ComfortScan,
Integrated hard drive and DVD-R/W,
Built-in memory for peripheral devices,
Imaging modes:
B-Mode,
Harmonic Imaging,
M-Mode,
Color Flow Mode (CFM),
Power Doppler Imaging (PDI) ,
PW Doppler,
Applications:
Obstetrician/Gynecologist,
Vascular,
Abdominal,
CEC,
Urology,
Cardiology,
Pediatrics,
Neonatal,
Transcranial,
Mobile on 4 wheels with brake,
Selectable menu in ENG, NOR, IT, FR, ESP, RUS , GREK, DUN, SWE, BR, NED, FIN, CHIN, JP,
Included:
User manual in English (PDF),
Power cord,
CONVEX GE-4C probe,
GE-12L linear probe,
GE-E8C endovaginal probe,
MITSUBISHI P93D Printer,
Dimensions: 52 x 62 x 136 cm,
Weight: 80 KG,
Cat 02,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

