



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



JUVENTAS gynecological chair JFG 2 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production,
Upholstery with medical certificate in dark brown color (upholstery can be replaced with a new one from our color chart for an additional fee),
Gynecological chair JFG 2 provides doctors and patients with full comfort during use. It was made of steel profiles, the structure of JFG 2 chair is solid and durable.
The set includes:
2x patient footrest,
Head rest,
Removable leg rest,
Patient bowl,
Features:
The chair has a fixed height,
Adjustable patient backrest,
Adjustable level of the entire chair,
Knee bolsters manually adjustable in all planes,
Removable leg rest,
Adjustable Trendelenburg and anti-Trendelenburg positions,
Trendelenburg tilt angle: from 0 to + 25 °,
Tilt angle of anti -Trendelenburg: from 0 to + 15 °,
Backrest tilt angle: from 0 to + 85 °,
Footrest adjustment in range: 38 cm,
Manual retractable step,
Dimensions:
Overall length: 186 cm,
Headrest: 36 x 24 cm,
Backrest: 66 x 90 cm,
Seat: 64 x 43 cm,
Leg rest: 65 x 48 cm,
Base: 55 x 105 cm,
Bowl: 30 x 15 x 7 cm,
Maximum load capacity 150 kg:
Weight: 86 kg,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

