



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



CHEST HI-801 Spirometer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Japanese production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Menu in language: French,
Specifications:
Flow rate:
Flow range: +/- 0.05 L/s. to +/- 14.0 L/s,
Flow resolution: 0.01 L/s,
Flow accuracy: +/- 3% L/s. or +/- 0.2 L/s.,
Volume:
Volume range: +/- 0.1 L. to 10.0 L,
Volume resolution: 0.01 L,
Volume accuracy: +/-3% or 0.05 L,
Screen:
Color display: 10.4 TFT,
Resolution: 640 x 480 px,
Built-in printer with 11.2 cm paper width,
Capable of storing up to 1000 patients on a memory card,
Includes:
Power cord,
SPIROMETER head,
Dimensions: 35 x 29 x 10 cm,
Weight: 4.5 KG,
Cat 01,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

