



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **CHEST HI-801 Spirometer (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Japanese production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Menu in language: French,  
Specifications:  
Flow rate:  
Flow range: +/- 0.05 L/s. to +/- 14.0 L/s,  
Flow resolution: 0.01 L/s,  
Flow accuracy: +/- 3% L/s. or +/- 0.2 L/s.,  
Volume:  
Volume range: +/- 0.1 L. to 10.0 L,  
Volume resolution: 0.01 L,  
Volume accuracy: +/-3% or 0.05 L,  
Screen:  
Color display: 10.4 TFT,  
Resolution: 640 x 480 px,  
Built-in printer with 11.2 cm paper width,  
Capable of storing up to 1000 patients on a memory card,  
Includes:  
Power cord,  
SPIROMETER head,  
Dimensions: 35 x 29 x 10 cm,  
Weight: 4.5 KG,  
Cat 01,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

