



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



PHILIPS PAGEWEWRITER TRIM II electrocardiograph ECG machine (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Manufactured by U.S.A.,
Power supply: 230 V,
Emergency power supply: through 1 built-in 12 V / 2.3 Ah battery,
Frequency: 50/60 Hz,
Software: 01.01,
Cardiograph
PageWriter TRIM II
by
Philips
is an advanced, reliable and compact ECG solution with ease of use,
Specifications:
12-lead ECG,
Alphanumeric keyboard,
Display: color TFT 640 x 480 px,
Sampling frequency: 2000 samples per second per electrode/duct,
Built-in memory for recording 50 tests,
Built-in PCMCIA card reader,
Report formats:
3x4 1R, 3R (Standard, Cabrera),
6x2 (Standard, Cabrera),
Panoramic 12 (Cabrera),
12x1 (Standard, Cabrera),
Rhythm (up to 12 selected leads),
12-lead extended measurements,
(1 minute of continuous waveform data for 1 selected lead),
Automatic frequency response:
0.05-150 Hz / 0.15-150 Hz / 0.5-100 Hz / 0.05-100 Hz / 0.15-100 Hz / 0.5-100 Hz / 0.05-40 Hz / 0.15-40 / 0.5-40 Hz,
Rhythmic frequency response:
0.05-150 Hz / 0.15-150 Hz / 0.05-150 Hz / 0.15-100 Hz / 0.05-40 Hz / 0.15-40 Hz,
Built-in thermal printer:
4-channel, real-time printing (paper width: 210 mm x 300 mm),
Menu in language: Dutch,
Included:
Power cord,
PHILIPS 10-lead ECG module,
ECG clips,
Dimensions:
Weight: 5 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

