



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDTRONIC LIFEPAK 20E defibrillator (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
U.S.A. production,
Power supply: 230 V,
Emergency power supply: via 1 built-in 11.1 V / 6 Ah battery
(no battery included)
- Can be retrofitted for an additional charge,
Frequency: 50/60 Hz,
Power: 120 VA,
Lifepak 20E
is a biphasic portable defibrillator, providing a complete response system for acute cardiac conditions,
Technical specifications:
AED mode,
Synchronization,
ECG parameters monitoring,
Easy-to-read 6.2" display,
Built-in printer (paper width: 50 mm),
Selectable menu in language: ENG,
Includes:
User manual in Polish (PDF),
Power cord,
ECG 3-lead cable,
AED cable,
Dimensions: 26 x 24 x 21 cm,
Weight: 6.5 KG,
Cat D133/2
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

