



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **ASHLEY RF micro-needle laser (Reconditioned)**

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Korean production, (2022),

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 150 VA,

The ASHLEY laser of the Korean company IDS, thanks to the use of isolated micro-needles, works thermally only at the set depth of the skin and does not damage the epidermis. This results in very good results in the treatment of stretch marks and scars,

The device is used for thermal, pinpoint damage to the dermis, with thermal damage to the epidermis omitted. This is possible, thanks to the tip with insulated or non-insulated needles automatically inserted into the skin to a set depth, administering RF current energy only below the skin surface, in a layer of about 0.5mm thick (depending on the set energy),

Features:

Ability to produce micro-damage to the skin at greater depths (up to 4mm), inaccessible to fractional lasers (up to 1-2mm),

Ability to produce skin microdamage at different depths in one treatment (e.g. 1mm + 1.5mm + 2mm),

Very short healing time and much lower risk of post-inflammatory hyperpigmentation,

Application:

Firming and improving skin tone,

Reduction of pores and wrinkles,

Scar reduction,

Specifications:

Two operating frequencies 1 and 2 MHz,

Depth of operation adjustable from 0.5mm to 4mm in 0.1mm increments,

Two modes of operation manual and automatic with single or double pulse,

Selectable menu in language: English,

Included:

Power cord,

Head,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for a period of 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

