



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



KAMPLEX AS7 diagnostic audiometer (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

U.K. production,

Power supply: 230 V,

Emergency power supply: through new built-in LR14 batteries,

Frequency: 50/60 Hz,

KAMPLEX AS7

is a portable audiometer that is ideal for schools, industry and primary medical practices, or wherever quick hearing screening is performed,

Specifications:

Frequency range: 250 - 8000 Hz,

Volume range: 0 to 80 dB,

Included:

Power cord,

Power supply,

AUDIOCUPS in-ear headphones,

Has up-to-date inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

