



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STIEGELMEYER PURO BREVO hospital bed (Ex-Display)

New (post-show),
Real photos of the product,
German manufacture, (2020),
Power supply: 230 V,
Emergency power supply: via built-in new battery,
Frequency: 50/60 Hz,
Hospital bed
PURO BREVO
With a low lying position and a large electric adjustment range of 28 - 86 cm,
PURO BREVO
convinces with its aesthetic appearance, excellent hygienic properties and effective protection against damage. The smooth surfaces exude minimalist elegance and high quality,
The compact design of the side rails gives the bed a slim appearance, integrated rollers and bumper cones will prevent damage to the bed and walls,
Technical specifications (electric adjustment with wired control remote control):
Trendelenburg tilt angle: 17°,
Anti-Trendelenburg tilt angle: 16°,
Backrest tilt angle: 70°,
Thigh support tilt angle: 40°,
Leg position elevation: 16°,
Bed height adjustment range: 28 - 84 cm,
Mobile bed with built-in central brake,
Includes:
Power cord,
Wired remote control,
Overall dimensions: 99.5 x 214.5 x 28/86 cm,
Maximum lifting capacity: 275 KG,
Weight: 155 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

