



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STIEGELMEYER PURO BREVO hospital bed (Ex-Display)

- New (ex-exhibition),
- **Illustrative photos,**
- German production, (2020),
- Power supply: 230 V,
- Emergency power supply: via a built-in new battery,
- Frequency: 50/60 Hz,
- **PURO BREVO** hospital bed with a low bed position and a wide range of electric adjustment in the range of 28 – 86 cm,
- **PURO BREVO** impresses with its aesthetic appearance, excellent hygiene properties, and effective protection against damage. Its smooth surfaces exude minimalist elegance and high quality.
- The compact design of the side rails gives the bed a slim look, integrated rollers and bumper cones prevent damage to the bed and walls,
- **Technical specifications (electric adjustment via wired remote control):**
 - Trendelenburg angle: 17°,
 - Anti-Trendelenburg angle: 16°,
 - Backrest angle: 70°,
 - Thigh rest angle: 40°,
 - Leg elevation position: 16°,
 - Bed height adjustment range: 28 – 84 cm,
- **Mobile bed with built-in central brake,**
- **Included in the set:**
 - Power cord,
 - Wired remote control,
- **Mattress not included - can be added for an additional fee,**
- Overall dimensions: 99.5 x 214.5 x 28/86 cm,
- Maximum load capacity: 275 KG,
- Weight: 155 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

