



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NDS LCD surgical monitor S.C.-SX19-A1511 (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Manufactured by U.S.A.,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 40 W,

The NDS LCD Surgical Monitor is an advanced high definition (HD) medical imaging system for minimally invasive surgery (MIS) and interventional procedures,

It incorporates NDS' proprietary Advanced Image Processing (AIP) technologies, providing a unique multimodal platform for simultaneous processing of images from different sources - endoscopic camera, fluoroscope, ultrasound or others - providing unparalleled accuracy, performance and collaboration for surgical and clinical teams,

Technical specifications:

Technology/Screen Type: a-Si TFT active matrix with in-plane switching,

Screen: 19" inches,

Display resolution: 1280x1024 (SXGA),

Viewing angle (horizontal/vertical): 178°,

Maximum brightness: 300 cd/m²,

Contrast ratio: 600:1,

Response time: 10-16 ms,

Connectors:

DVI-1 (IN/OUT) / VGA IN / S-VIDEO (IN/OUT) / HD SDI (IN / OUT) / COMP-SOG / DVI-2 IN RGBS / SERIAL CONTROL / R-G-B / HS C-SYNC / VS / GPIO, Y / C,

Selectable menus in ENG, GER, FR, IT, SV, ESP, NED,

Included:

Power cord,

Dimensions: 18 x 4 x 16 cm,

Weight: 6.8 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

