



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ARIATELE SPACELABS ECG telemetry transmitter (Reconditioned)

Reconditioned (used),
Technical condition: good,
Visual condition: very good,
Real photos of the product,
Manufactured by U.S.A., (2019),
Powered by 2 new built-in AA batteries,
Frequency: 50/60 Hz,
Telemetry transmitter
ARIATELE SPACELABS ECG
Has an easy-to-read color graphic display. Up to four ECGs can be viewed along with heart rate,
Features:
Basic ECG and multi-parameter LCD display,
Color display,
Detection of pacemaker heart rate display,
Specifications:
Display: color LCD,
The pacemaker detection function detects pacemaker pulses of ± 0.5 mV with a pulse width of 0.2 to 2 ms and ± 1.0 mV with a pulse width of 0.1 ms,
Pacemaker rejection pulse period of 0.1 ms width, amplitude up to 700 mV (tested at 2 mV, 10 mV, 300 mV and 700 mV), with an overshoot of no more than 4 ms,
Included:
5-lead ECG lead,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

