



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MAQUET operating table (Reconditioned)**

Product temporarily unavailable

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product

, German manufacture,

New upholstery

with medical certificate in black

, 5 segmented,

MAQUET operating table

is an easy-to-use, manual operating table designed for many surgical fields,

Adjustable manually in many planes:

Height in range: 80 - 100 cm,

Backrest angle in the range of -7.5° to +60°,

Lateral tilt: ±25°,

Trendelenburg position: 20°,

Anti-Trendelenburg position: 20°,

Extendable headrest: + 15 cm,

Features:

X-Ray translucency,

Rails for attaching additional accessories,

Stainless steel base,

Mobile with brake,

Possibility to retrofit the table (for an extra charge) with additional accessories, examples shown in the last photos,

Dimensions:

Bed: 188 x 47 cm,

Base: 75 x 97 cm,

Weight: 220 kg,

Has a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

3 months for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan,

