



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## ALCON phacoemulsification handpiece set (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- **Included in the set:**
  - Handle: **INFINITI OZIL 8065750469,**
    - **Clinical use:**
      - Emulsification of hard cataract nuclei,
      - Reduction of thermal burns at the cutting site,
  - Handle: **ACTIVEENTRY AS 8065752914 (3 pcs.):**
    - **Clinical use:**
      - High-risk cataract surgery (including hard nuclei, weakened ligaments),
      - Procedures requiring precise control of chamber depth,
  - Handle: **INTREPID AUTOSERT 8065751755:**
    - **Clinical use:**
      - IOL implantation through micro incisions ( $\leq 2.2$  mm.),
      - Minimally Invasive Surgery (MICS),
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
  - **6 months for domestic market (Poland),**
  - **3 months to international market,**
  - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

