



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS TC-C1 endoscope cart (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
U.K. manufacture,
OLYMPUS TC-C1
is a compact endoscope cart for transporting endoscopes and delicate medical equipment,
Equipped with:
Transport handle,
Rubber bumpers,
3 shelves with dimensions: 39 x 45 cm,
Lower top with dimensions: 37.5 x 45 cm,
Camera head holder,
Made of powder coated steel,
Mobile on 4 casters with a diameter of 10 cm including 2 with a brake,
Dimensions: 117 x 56.5 x 48.5 cm,
Weight: 32 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

