



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## MELAG TYPE 28 Autoclave (Reconditioned)

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 1500 W,  
MELAG TYPE 28  
Is a steam autoclave,  
Features:  
Built-in thermometer and pressure gauge,  
Indicator lights for autoclave start-up and heating,  
German language control panel,  
Water infusion and quick pressure release function,  
Possibility to set a timer for max. 30 min,  
Water drainage connector on the front of the device,  
On the lid of the autoclave a place to put the tray,  
All steam guiding elements made of stainless materials,  
Specifications:  
Capacity of sterilization chamber: 14.51 L,  
Pressure range: 2 bar (134o C) - 2.5 bar (138o C),  
Maximum load capacity: 2.5 kg of instruments including trays or 0.2 kg of textiles,  
Water consumption during sterilization: 400 ml,  
Has 1 sterilization program: Pressure: 2.5 Bar - time: 30 min - temp: 134° C max,  
Includes:  
2x sieve tray with dimensions: 28 x 18 cm,  
Tray holder,  
Stand with removable handles for 4 trays,  
Instruction manual in German (PDF),  
Dimensions:  
General: 57 x 44.5 x 46 cm,  
Compartments: Ø 23 cm, depth 41 cm,  
Weight: 20 KG,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

