



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **MGE SAM 12 surgical suction machine (Ex-Display)**

New (post-show),  
Real photos of the product,  
U.K. production, (2020),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 90 W,  
Surgical suction machine  
MGE SAM 12

is an ideal choice for both hospitals and private health centers. It is a suction unit providing high vacuum and a maximum flow of up to 30 l/min,

The suction unit is equipped with a modern, maintenance-free and also oil-free pump, thanks to which it provides low noise levels and does not emit harmful oil vapors. Construction made of aluminum alloy guarantees low weight, making the device lightweight and comfortable to transport,

Specifications:

Maximum flow rate up to: 30 l/min,

Maximum pressure:

0 - 760 mm/Hg,

0 - 0.85 bar,

Loudness of operation: 50 dB,

Included:

Power cord,

Autoclavable jar (121° C) SERRES 2 l. capacity with disposable cartridge,

Hose for the jar,

Biological filter,

Overall dimensions: 45 x 20 x 38 cm,

Weight: 10 KG,

It has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

