



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **JULABO SW22 shaking water bath (Reconditioned)**

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

German production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 2200 W,

Water baths

SW22 JULABO

are modern devices ideal for everyday use in the laboratory. Their robust design and stainless steel construction provide them with extreme strength and durability, as well as stability even at the highest shaking speeds,

Features:

Sloped tank edge,

Liquid level control,

Audible alarm signal in case of failure,

Removable shaking platform,

Spill-proof and splash-proof design,

Simple programming,

Specifications:

Temperature range 20 to 99.9 °C,

Shaking amplitude: 15 mm,

Rotational speed: From 20 to 200 rpm,

Temperature stability: +/- 0.2,

Included:

ENG user manual (PDF),

Power cord,

Dimensions:

Inner bath: 50 x 30 x 18 cm,

External: 70 x 35 x 26 cm,

Total weight: 20 KG,

It has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

