



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEMMERT U15 laboratory warmer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 1400 W,
MEMMERT U15
is a capacious and easy-to-use laboratory hothouse,
Equipped with:
Indicator lights indicating that the device is started and the heater is working,
Knob for temperature control in the range: 50° - 220° C,
Built-in thermometer with a range of 0° - 250° C,
Knob with time control,
Two perforated trays with adjustable height of their position,
6-stage ventilation control,
Includes:
Power cord,
2x trays included,
Dimensions:
Exterior: 56 x 42.5 x 59 cm,
Compartments: 46 x 25.5 x 26 cm,
Trays: 22.5 x 47 cm,
Weight: 22 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

