



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **EMS MEDI-LINK 70 short-wave diathermy with ultrasonic module 72 (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
U.K. production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 200 VA,  
Software: V1.12,  
Device  
EMS MEDI-LINK MODEL 70  
is a multifunctional physical therapy device, mainly used in rehabilitation, pain management and treatment of musculoskeletal disorders,  
Equipped with a module  
MEDI-LINK 72:  
Frequency: 1.1 / 3.3 MHz,  
Maximum output power: 10 W,  
Maximum intensity: 2W/cm  
2  
,  
Adjustable treatment time in the range of 30 minutes,  
CW and pulsed modes 1:1 / 1:2 / 1:4 / 1:9,  
Application:  
Deep tissue penetration - accelerates healing, reduces inflammation and swelling,  
Analgesic effect - arthritis, tendonitis,  
Scar softening - helpful after surgeries or injuries,  
Selectable menu in English, French,  
Included:  
User manual in English (PDF),  
Power cord,  
EMS head 0.5 cm  
2  
,  
EMS head 5 cm  
2  
,  
Dimensions: 33 x 20 x 10 cm,  
Weight: 3 KG,  
It has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

