



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Infusion (syringe) pump B.BRAUN PERFUSOR FM (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Emergency power supply: via built-in new battery,
Frequency: 50/60 Hz,
Power: 18 VA,
Syringe selection: 10/12, 20, 50/60 ml,
Automatic recognition of syringe size,
Dosage: 0.1 - 200 ml/h,
Preselection dose: 1 - 999.9 ml,
Bolus function settings:
Rate: 1 ml/h to 1800 ml/h,
Volume: 0.1 ml to 50.0 ml,
Interval time: 1 min to 99 h 59 min,
Alarms:
Battery alarm / battery pre-alarm,
Syringe alarm / syringe pre-alarm,
Pressure alarm,
Standby alarm,
Has two displays, liquid crystal and LED,
FM (fluid manager system) socket,
Equipped with a carrying handle,
Menu in language: German,
Includes:
User manual in Polish (PDF),
Power cord,
Dimensions: 33.5 x 19 x 11 cm,
Weight: 3.5 KG,
Cat 02,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

