



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **BIOMEDIC COMFORT LINE dialysis chair (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 173 VA,  
Material upholstery with medical approval in blue color (upholstery can be replaced with a new one from our color chart for an additional fee),  
Dialysis chair  
**BIONIC COMFORT LINE**  
provides free adjustment with wide-range positioning,  
The chair has excellent mobility making it suitable for free movement around medical facilities  
Technical specifications:  
Adjustable by wired handheld remote control:  
Seat height adjustment in the range: 60 - 80 cm,  
Backrest adjustment in the range: 75°,  
Leg rest adjustment within: -30° / 0° / 5°,  
Trendelenburg tilt in the range of 22°,  
Electrically extendable patient step within 30 cm,  
Programmed button with the ability to adjust the chair to 0 "zero" position,  
Mobile chair with built-in handbrake on both sides,  
Freely adjustable patient armrests on both sides,  
Includes:  
Wired remote control,  
Dimensions:  
Headrest: 44 x 36 cm,  
Backrest: 62 x 64 cm,  
Seat: 56 x 62 cm,  
Leg rests: 52 x 62 cm  
Overall length: 190 cm  
Maximum load capacity: 240 KG,  
Weight: 90 KG,  
Cat 02,  
Has a current inspection and is ready for use,  
Issued Technical Passport (Service Report) valid for a period of 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

