



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BIOMEDIC COMFORT LINE treatment and cosmetic chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 173 VA,
Material upholstery with medical certificate in blue color (upholstery can be replaced with a new one, from our color pattern for an additional fee),
Treatment and cosmetic chair
BIONIC COMFORT LINE
provides free adjustment with the possibility of positioning in a wide range,
The chair has excellent mobility due to which it is adapted for free movement around medical facilities
Technical specifications:
Adjustment by wired handheld remote control:
Seat height adjustment in the range: 60 - 80 cm,
Backrest adjustment in the range: 75°,
Leg rest adjustment within: -30° / 0° / 5°,
Trendelenburg tilt in the range of 22°,
Electrically extendable patient step within 30 cm,
Programmed button with the ability to adjust the chair to 0 "zero" position,
Mobile chair with built-in handbrake on both sides,
Bilateral / freely adjustable patient armrests,
Includes:
Wired remote control,
Dimensions:
Headrest: 44 x 36 cm,
Backrest: 62 x 64 cm,
Seat: 56 x 62 cm,
Leg rests: 52 x 62 cm
Overall length: 190 cm
Maximum load capacity: 240 KG,
Weight: 90 KG,
Cat 02,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

