



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS WM-NP1 endoscope cart (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
U.K. production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Olympus WM-NP1
is a stable and sturdy cart suitable for carrying endoscopes and delicate medical equipment,
It has:
Upper top with dimensions: 47 x 38 cm,
Lower top with dimensions: 47.5 x 56 cm,
2 shelves with dimensions: 40 x 54 cm with manual height adjustment (every 1.5 - 2 cm) with Allen screws,
Handles for easy movement of the cart,
Holder for endoscopes with manual height adjustment of 86 cm,
Electrical installation, 12 power sockets for transported medical equipment,
Product features:
Power switch on the front of the cart,
Detachable rear panel of the cart for easy cable management,
Power cords aesthetically concealed in panels on the sides of the cart,
Corners of shelves protected with protective caps and edges over wheels with rubber,
Mobile on four wheels with Ø 12 cm, two of which have brakes,
Included:
Power cord,
3 pcs power cord for the device,
Wiring:
VGA monitor,
SDTV: SYNC/B/R/G/COMP,
HDTV: SYNC/B/R/G,
MONITOR SVHS IN,
VCR SVHS:
IN,
OUT,
VCR RS-232,
DIGITAL CAPTURED SVHS IN,
VIDEO PRINTER OUT,
PHOTO: SYNC/B/R/G,
Dimensions:
General: 67 x 66 x 133 cm,
Weight: 87 kg,
Cat 08,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

