



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



WASSENBURG MEDICAL endoscope washer WD440 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production, (2017),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 3 kVA,
Washer
WASSENBURG WD440

are designed for a controlled, fully automatic process of cleaning and disinfection of flexible endoscopes and other thermolabile devices. Their design and functionality ensure the highest standards of safety and efficiency in an unprecedentedly short time (a full process in as little as 22 minutes),

WASSENBURG WD440 has a unique contoured chamber with a corrugated bottom connected to a special spiral tube. As the solution enters the chamber and flows through the tube, it gently lifts the endoscope along with its probe, eliminating contact points and allowing efficient penetration of any area. The basis of the system's effectiveness is proven pump and connection technology. Each endoscope channel is independently processed and controlled for flow, pressure and leakage,

Two independently operating separate washing chambers provide easy access to individually loaded and connected endoscopes and increase the efficiency of the washer. The chambers are made of the highest quality stainless steel for optimum performance of the cleaning and disinfection process over the life of the machine. In addition, their corrugated bottom minimizes contact with the endoscope. The chambers are designed to reduce water consumption as much as possible and ensure that the process temperature remains constant,

Specifications:

Water consumption: 35 l. of water per 1 standard cycle,

Screen: 7-inch (TFT touch screen),

Min / Max Pressure: 2.5 / 4 bar,

Min. water flow: 10 l/min,

Noise level: 60 - 65 dBA during use, 50 dBA at rest,

Process time: 22 minutes,

Included:

Power cord,

GRYPHON wired reader,

Can be retrofitted with connections for the endoscope of your choice (at additional cost),

Overall dimensions: 94 x 75 x 103 cm,

Weight: 200 KG,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

