



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybierański
+48 882 762 006

kw@praiston.pl



Surgical diode laser SWING 15 1470 nm + tumescent anesthesia pump NOUVAG AG DISPENSER DP30 (Ex-Display)

New (post-show),
Real photos of the product,
Polish production, (2021),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 350 VA,
Laser
SWING
Emits radiation with a wavelength of 1470 nm. This is the optimal wavelength for penetration into the skin. The correspondingly low absorption in hemoglobin and melanin, while high for selective photothermolysis in water and fat,
These features make the laser
SWING 15
excellent tools in endovascular, ENT, laser lipolysis, as well as a number of other procedures. Importantly, the radiation spectrum is safe for vision.
Features:
Painless procedures,
No risk of damage to adjacent tissues,
Minimal risk of complications,
Outpatient procedure as part of same-day surgery,
Quick return of the patient to activity,
No need for incision,
Procedure without the use of rubber bands, clamps and threads,
No need for suturing,
Minimization of bleeding,
Possibility of repeating the procedure,
Eye-safe radiation spectrum,
Intuitive interface,
Large, convenient touch screen,
Application:
Aesthetic medicine
- Lipolysis, laser-assisted liposuction,
Phlebology
a - endovascular removal of insufficient venous trunks (EVLT),
Proctology
- Ablation of grade II, III and IV hemorrhoids,
Laryngology
- Microsurgery of tumors, polyps and strictures in the larynx, laser palate and uvula plasty, treatment of snoring, excision of tonsils and adenoid growths, removal of nasal polyps and adhesions, reduction of hypertrophied nasal auricles, oral surgery, outpatient myringotomy, middle ear surgery,
Gynecology
- Hysteroscopy and laparoscopy, resection of endometriosis foci, bladder excision, myomectomy, polypectomy, open surgery (coagulation of hemangiomas),
Dermatology,
Specifications:
Output power: 15 W,
Wavelength: 1470 nm,
Screen: color (touchscreen),
Fiber optic diameter: 400 µm / 600 µm,
o
Laser equipped with
NOUVAG AG DISPENSER DP 30 tumescent anesthesia pump
:
DISPENSER DP 30
is a universal pump used for tumescent anesthesia in endovascular, liposuction and other procedures,
Intuitive to use,
Equipped with a metal stand for hanging a bottle of tumescent fluid,
Infinitely precise flow rate control with potentiometer from 0 - 100% ,
Foot switch allows you to turn the flow on and off and control the discharge speed without interrupting the operation,
Specifications:
Flow rate: 1 - 15 L/h,
Maximum pressure: 2 bar,
Includes:
NOUVAG Model VARIO foot attachment,
Dimensions:
Laser: 43 x 46 x 20 cm,
Pump: 25.5 x 23 x 10 cm,
Weight
:
Laser: 13 KG,
Pumps: 3.5 KG,
It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

