



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



BTL EXELIS ULTRA 360 tissue contouring system (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
U.K. production, (2022),
Power supply: 230 V,
Frequency: 50/60 Hz,
BTL Exilis Ultra 360
is a treatment that combines two technologies - ultrasound with high-frequency electromagnetic wave (radiofrequency, RF),
It involves the controlled heating of adipose tissue to a temperature of 40 °C, resulting in both improved formation of new collagen fibers and fat cell reduction,
This device models and firms the body "from head to toe". The advanced system for non-invasive aesthetic treatments can be applied to as many as 20 areas of the face and body,
Application:
Fat reduction (pronounced fat loss in treated areas),
Fat volume reduction,
Improving body contour
Silhouette slimming,
Remodeling of selected treated body areas (especially areas resistant to diet and exercise),
Firming the skin and improving skin tone (especially after pregnancy),
Increasing skin density,
Reduction of cellulite,
Reduction of stretch marks,
Included:
Power cord,
R1 head,
R2 head,
Weight: 30 KG,
Cat 02,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

