



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAQUET RADIUS gynecological chair 1557.01AC (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 290 VA,
Maquet Radius
is used for gynecological and urological examinations and procedures,
Upholstery will be replaced with a new one in a color chosen by the customer from our color chart,
Possibility of electromechanical adjustment by foot control and wired handheld remote control:
Back segment angle in the range: 57°,
Seat height adjustment in the range: 61 - 91 cm,
Features:
Storage compartment for a 38 cm wide roll of backing paper,
Cutout in the upholstery of the back segment for backing paper,
3 pre-programmed buttons on handheld remote control,
S button to stop the programmed button mode,
Pull-out waste bowl with Ø 21 cm,
Ability to electrically unfold the chair into the reclining position,
Includes:
User manual in Polish (PDF),
Power cord,
Wired foot remote control,
Wired handheld remote control,
Waste bowl,
Dimensions:
Back segment: 91 x 61 cm,
Seat: 61 x 31 cm,
Length in extended position: 175 cm,
Overall: 155 x 76 x 138/168 cm,
Maximum load capacity: 180 KG,
Weight: 94 KG,
Has up-to-date inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

