



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CACAN 800 ophthalmological chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Dutch production, (2013),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 650 W,
Upholstery with medical certificate in blue color (possibility to replace the upholstery with a new one from our color pattern for an additional fee),
CANAN 800 ophthalmic chair
is a multifunctional chair ideally suited for ophthalmology and ENT surgeries, and through built-in functions can be used as a treatment chair,
Technical specifications:
Electric control of the chair through built-in foot joystick on both sides:
Backrest tilt in the range: 90° / 0° / -3° (automatic synchronization of leg rest),
Seat height adjustment in the range: 66 - 97 cm,
Programmed electric adjustment of selected positions in 5 programs,
Possibility to electrically unfold the chair in the reclining position,
Patient armrest on both sides,
Free tilt adjustment of the headrest,
Dimensions:
Backrest: 44 x 52 cm,
Seat: 48 x 44 cm,
Leg rest: 44 x 54 cm,
Total length in lying position: 168 cm,
Weight: 125 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

