



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **ENT chair CACAN 800 (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Dutch production, (2013),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 650 W,

Upholstery with medical certificate in blue color (it is possible to replace the upholstery with a new one, from our color pattern for an additional fee),

ENT chair CACAN  
800

is a multifunctional chair ideally suited for ophthalmology and ENT surgeries, and through built-in functions can be used as a treatment chair,

Technical specifications:

Electric control of the chair through built-in foot joystick on both sides:

Backrest tilt in the range: 90° / 0° / -3° (automatic leg rest synchronization),

Seat height adjustment in the range: 66 - 97 cm,

Programmed electric adjustment of selected positions in 5 programs,

Possibility to electrically unfold the chair in the reclining position,

Patient armrest on both sides,

Free tilt adjustment of the headrest,

Dimensions:

Backrest: 44 x 52 cm,

Seat: 48 x 44 cm,

Leg rest: 44 x 54 cm,

Total length in lying position: 168 cm,

Weight: 125 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

