



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FAMED ŻYWIEC SZ-01 operating table (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production,
Famed SZ-01
is a table whose simple and robust design allows it to be used also outside the operating room,
Upholstery with medical certificate in black color,
Equipped with:
Head segment mounted on a rail located above the back segment,
Removable, split leg segment to expand the range of treatment options,
Rails for attaching additional accessories, including one above the head segment,
Locking angle adjustment of the back segment,
Two comfortable handles on both sides of the head segment for convenient table guidance, including one for setting the table in Trendelenburg and anti-Trendelenburg positions,
Hydraulic control of the height of the bed with a foot lever in the range of 81 - 106 cm,
Infinite manual adjustment of the angle of the sleep surface and its individual segments by means of a gas spring lever:
Change of the patient's lying direction relative to the table column: Trendelenburg position in the range of 25° and anti-Trendelenburg position in the range of 12°,
Segment:
Back : (-40)° / 0° / 80°,
Head : (-40)° / 0° / 50°,
Legs (split): (-90)° / 0° / 15°,
Bed height adjustment within: 23 cm,
Ability to visit both elements of the split leg segment in the range of 0 - 90°,
Transparent for X - RAY radiation,
Adapted to C-arm,
Mobile, on Ø 12.5 cm wheels with central brake,
Dimensions:
Bed: 55 x 202 x 5 cm,
Overall: 200 x 56 x 81/106 cm,
Maximum load capacity: 160 kg,
Weight: 117 kg,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

