



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAQUET ALPHASTAR operating table 1132.12A2 (None)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in batteries 12 V / 17 Ah,
- Frequency: 50/60 Hz,
- Power: 400 VA,
- **New upholstery (memory visco) with medical certificate in black,**
- **Adjustable in many planes using a wired remote control or via the panel on the column:**
 - ON/OFF,
 - LOCK/UNLOCK - locking and unlocking the table,
 - NORMAL/REV - patient orientation on the operating table,
 - TREND/REV - change of the patient's lying direction relative to the table column: Trendelenburg and anti-Trendelenburg position within the range of $\pm 30^\circ$,
 - The table is equipped with an auxiliary drive controlled by a handheld remote control, allowing forward and reverse travel,
 - UP/DOWN - height of the operating surface (up/down),
 - TILT L/TILT R - lateral tilt in the range of 20° ,
 - BACK UP/BACK DN - adjustment of the back section position,
 - LEG UP/LEG DN - adjusting the leg segment: up/down,
 - FLEX/REFLEX - bending the tabletop in half in both directions: roof position / V position (Beach chair),
 - LEVEL - reset: return to the horizontal position of the bed,
- **Electric adjustment via wired foot control:**
 - TREND/REV - change of the patient's lying direction relative to the table column: Trendelenburg and anti-Trendelenburg position within the range of $\pm 30^\circ$,
 - TILT L/TILT R - lateral tilt in the range of 20° ,
 - UP/DOWN - height of the operating surface in the range: 68 – 112 cm (up/down),
- **Manually adjustable:**
 - Headrest,
 - Leg rest segment,
- The base housing is made of durable plastic,
- The table has X-RAY transparency,
- Equipped with rails for mounting additional equipment,
- Table adapted to C-arm,
- Mobile, on double wheels
- Manually removable headrest and foot section,
- **Included in the set:**
 - Wired remote control,
 - Wired foot controller,
 - Power cord,
- **The table can be equipped with other accessories from our offer for an additional fee, examples are shown in the last photos,**
- **Dimensions:**
 - Headrest: 54 x 26 cm,
 - Back support: 54 x 30 cm,
 - Bed segment: 54 x 100 cm,
 - 2-segment leg rests: 25 x 80 cm,
 - Total length: 242 cm,
- Maximum load capacity: 450 KG,
- Weight: 325 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

