



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Surgical diathermy (argon) ERBE ICC-200 with APC-300 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 200 W,
Technical specifications:
Infinitely adjustable power,
Digital displays of set power,
Foiled, touch-sensitive control panel,
Features the following functions:
ENDO-CUT,
Autocut - max. 200 W,
Autocoagulation - max. 120 W,
Available coagulation modes: SOFT, FORCED, BIPOLAR,
Included:
ERBE APC 300 argon coagulation attachment:
Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German manufacture,
Power supply: 230 V,
Frequency: 50/60 Hz,
Adjustable gas (argon) flow from 0.1 l/min to 9 l/min
Eight work programs,
Large electronic easy-to-read display,
ERBE mobile diathermy cart:
Storage compartment for two 5-liter argon cylinders,
Four swivel casters (2 with brake),
Includes:
Foot attachment,
Factory new bipolar cable, BOWA forceps (bipolar tweezers can be retrofitted at additional cost),
Factory new monopolar electrode lead with BOWA grip - two buttons, 4 mm shaft (no electrodes included, can be retrofitted at additional cost),
Factory new cable with BOWA clip for disposable neutral electrodes (no electrodes included, can be retrofitted at additional cost).
,
User manual in English (PDF),
Power cord,
Pressure regulator with pressure gauge and gas hoses,
2x argon cylinder (5 l.),
Overall dimensions: 45 x 47 x 120 cm,
Weight: 75 KG,
It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

