



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PRIM medical chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production,
Upholstery with medical certificate in yellow (possibility to replace the upholstery with a new one from our color pattern),
Height adjustable with a gas actuator,
Infinitely adjustable seat height in the range: 51 - 71 cm,
Chrome base with diameter: 50 cm,
Mobile on double casters,
Ability to rotate the seat around its own axis in the range of 360°,
Dimensions:
Seat diameter \varnothing : 41 cm,
Weight: 20 KG,
It has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

