



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **TOF-WATCH neuromuscular conduction monitor (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Irish production,  
Powered by 1 built-in new 9-volt battery,  
TOF-Watch  
is a microprocessor-controlled, easy-to-use , portable neuromuscular conduction monitor using the principles of acceleromyography,  
It provides precise information about the degree of flaccidity,  
Enables faster recovery of normal muscle and respiratory function,  
Prevents the formation of persistent neuromuscular block,  
Automatically shuts off after 2 h of inactivity,  
Features an easy-to-read 3" inch display,  
Stimulation modes:  
Quadruple series (TOF),  
Power number (PTC),  
Single pulse (ST) 1 or 0.1 Hz,  
Double pulse (DBS) 3.3 or 3.2,  
Tetanic stimulation (TET) 50 or 100 Hz,  
Specifications:  
Stimulation current: 0 - 60 mA at impedance less 5 kOhm,  
Single-phase pulse duration: 200  $\mu$ s,  
Included:  
3-lead patient cable,  
Dimensions: 20 x 9.5 x 5 cm,  
Weight: 1 KG,  
Cat 03,  
Has a current inspection and is ready for use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

