



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



WHALEDENT UC1-230 ultrasonic cleaner (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Manufactured by U.S.A.,

Power supply: 230 V,

Frequency: 50/60 Hz,

Ultrasonic cleaning is based on the generation of high-frequency sound waves. The sound waves, aided by specially formulated cleaning solutions, produce a high-energy cavitation effect. During cavitation, millions of tiny bubbles are formed, which then collapse (implode), releasing huge amounts of energy and shock waves, which in turn scrub the surface of dental instruments, equipment and other apparatus placed in the cleaning solution. The powerful scrubbing action reaches the tiniest crevices for which manual brushing would be unattainable. The combination of energy with specially formulated solutions makes ultrasonic cleaning the most effective method for removing both large and microscopic waste,

Specifications:

Simple and intuitive operation,

Housing made of high quality plastic,

Working frequency: 55 kHz,

Adjustable time from 0 to 15 minutes and continuous operation,

Includes:

Basket,

Dimensions:

General: 34 x 25 x 23 cm,

Compartments: 29 x 15 x 10 cm,

Basket: 25.5 x 10.5 x 9.5 cm,

Weight: 3.5 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

