



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **GOWLLANDS MAYFIELD Ophthalmoscope (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
U.K. production,  
Power supply: 230 V,  
Light source powered by 2x LR14 battery,  
Frequency: 50/60 Hz,  
Ergonomic shape of the sight head, easily accessible apertures, filters and dioptic dial allow changing parameters without interrupting the examination,  
The LED illumination is finely tuned to the optical system and as a result provides exceptionally bright and even light for easy diagnosis,  
Specifications:  
Diopter adjustment range: + 25D to - 40D,  
Included:  
ENT mirror,  
Nasal speculum,  
It has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,  
The warranty period does not include bulbs and other light sources installed in the device,  
Possibility of financing (Poland only): Installments, Leasing, Loan,

