



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **CPM rehabilitation splint ORMED ARTROMOT K2 PRO-CHIP (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 40 VA,

Rail  
CPM Artromot K2  
PRO-CHIP

is a modern and advanced medical device designed for rehabilitation of the knee and hip joint. Dedicated to patients after surgeries, injuries and orthopedic procedures, it is the ideal solution to support a quick and safe return to full mobility,

### **Features:**

#### **Rehabilitation Effectiveness:**

The Artromot K2 PRO-CHIP CPM splint is designed to provide painless, automatic and controlled mobilization of the knee and hip joint. The CPM (Continuous Passive Motion) system allows smooth, uninterrupted joint movement without patient effort, which stimulates tissue regeneration, improves circulation and reduces swelling,

#### **Ease of use:**

The device is intuitive to use with a simple control panel. Any user, whether patient or medical staff, will quickly learn how to use the Artromot K2, making daily rehabilitation significantly easier,

#### **Individual adjustment:**

The Artromot K2 PRO-CHIP splint is suitable for patients of different sizes and degrees of conditions. The device allows for precise adjustment of the joint's range of motion, allowing for gradual and safe rehabilitation, tailored to the individual patient's needs,

#### **Advanced technology:**

The Artromot K2 PRO-CHIP CPM splint has state-of-the-art features that allow monitoring the progress of therapy and adjusting parameters according to rehabilitation needs. As a result, the device supports comprehensive rehabilitation therapy, ensuring rapid recovery,

#### **Safety:**

The splint is made of high-quality materials that ensure its strength and durability, even with long-term use. Artromot K2 also has numerous certificates of compliance with European medical safety standards,

#### **Application:**

Post-operative rehabilitation after arthroscopy, ligament reconstruction, knee and hip endoprosthesis,

Treatment of post-traumatic conditions such as fractures or dislocations of the knee and hip joints,

Prevention of adhesions and contractures that can occur after orthopedic surgeries,

Rapid return to mobility by stimulating the body's natural regenerative processes,

#### **Specifications:**

Range of motion in the knee joint: -5° / 0° / 115°,

Range of motion at the hip joint: 0° / 10° / 110°,

Free speed intensity adjustment in percentage range from 1% to 100% ,

Adjustable speed of movement and the ability to set pauses at selected points of movement (pause range up to 30 seconds),

Adjustment of treatment execution time in the range: 1 to 300 minutes and continuous operation,

Function of automatic joint warm-up at the beginning of therapy, which increases the effectiveness of rehabilitation,

Mobility - lightweight and compact design allows for easy transportation and use both in hospitals, clinics and patients' homes,

The device is designed for patients in the height range of 135 - 195 cm,

#### **Includes:**

Power cord,

Wired remote control,

Overall dimensions: 32 x 56 x 100 cm,

Weight: 14 KG,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

#### **Warranty:**

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

