



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



TOPCON RM-8000B autorefractometer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Japanese production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 60 VA,
RM-8000B
is a device for objective refraction measurement among Topcon autorefractometers,
It has the following diagnostic capabilities:
Refraction:
Contact lens measurement option,
Patient PD,
Patient's VD - 0/12/13.5 mm,
Cylinders (CYL): 0 to $\pm 8D$ in 0.25D step (0.12D step available),
Sphere (SPH): -25 to + 22D in 0.25D step (0.12D step available),
Axis: 1 ° to 180 ° in 1 ° step (5 ° step available),
Minimum lens diameter: 2mm,
Fixation: Autofocus system,
Features:
Built-in printer,
Manual canopy adjustment,
5 inch screen,
Outputs for computer connection,
Menu in language: English,
Included:
Power cord,
Overall dimensions: 27 x 48 x 45 cm,
Weight: 19 KG,
Has up-to-date inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

