



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS OTV S7PRO HD 12E camera head (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Japanese production,

Olympus OTV - S7 PRO

is a camera head that provides high resolution and extremely accurate imaging with its 1080/50i PAL format,

Narrow band imaging

(NBI)

is an optical image enhancement technology that improves the visibility of vessels and other tissues on the mucosal surface. The technology takes advantage of differences in the absorption and reflection properties of blood and tissues surrounding a blood vessel, making it easier to image capillaries that contrast with the surrounding mucosa,

Technical specifications:

Magnification 1.2x,

The image can be magnified up to 1.5 x from the original size while maintaining high HD image quality,

It has 3 switches built to program frequently used functions, such as electronic zoom and image capture, for convenient remote control,

Sealed camera head housing prevents lens fogging after steam sterilization,

Sensor system: Interline type solid-state CCD image sensor,

Focal distance (f): 15.8 - 31.3 mm,

Included:

OLYMPUS container for sterilization,

Dimensions:

Camera head: 11 x 6.5 x 4 cm,

Cable length: 4 m,

Sterilization container: 48 x 23 x 9 cm,

Weight: 0.5 KG

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

