



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SILENT AIR DA7001 30L Dental Compressor (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Manufactured by P.R.C, (2017),

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 750 VA,

Features:

Oil-free and ultra-quiet design provides clean and harmless air supply. There is an anti-corrosion coating inside the tank,

Copper motor provides durability, low temperature and low noise,

With electric protection device and rust-proof management, safe operation can be ensured, and the single motor can be well controlled for regular operation,

The motor is made entirely of copper, which ensures 24-hour continuous operation of the whole machine,

The anti-corrosion and sterilization coating in the gas tank ensures the purity of the gas source,

Advanced multiphase filtration ensures exceptionally clean air and prolongs service life,

Original oil-free piston with high efficiency and low dispersion that produces clean air to meet the needs of connected equipment,

Cylinder sleeve that has been hard anodized to increase its life by 40% compared to conventional sleeves,

Hard and sturdy housing to ensure long service life of the device,

Specifications:

Power: 0.75 kW,

Maximum air flow: 152L/min,

Air flow at 7 bar: 78L/min,

Tank: 30L,

Noise ≤ 68 dB(A),

Dimensions: 41 x 41 x 65 cm,

Weight: 27.6 KG,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

