



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **ATTOS diode laser (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
P.R.C. production, (2021),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
ATTOS

is a diode laser for hair removal designed for small beauty salons that value economical solutions,  
The head is equipped with a treatment spot of 12 mm x 20 mm, which allows to effectively depilate even large parts of the body such as legs, back or chest. The output power of the laser guarantees the effectiveness of the epilation treatment. A complete treatment includes 3 to 6 treatment sessions,

The ATTOS diode laser is equipped with an efficient TEC Cooling system. The device uses up to four cooling methods: sapphire glass, water, air and semiconductor. The system cools the head to 21 degrees Celsius (frost appears on the head while the device is operating). This effectively lowers the skin temperature and makes the treatment comfortable,

Wavelengths:

1064 nm - ideal for removing light hair from dark skin,

808 nm - the beam with the most effective effect on melanin, best for black and dark hair,

755 nm - extremely effective on fine, thin hair with low melanin content, guaranteeing low soreness of the treatment,

Specifications:

Energy: 10 - 150 J/cm

2

,

Pulse length: 10 - 800ms fully adjustable,

Frequency: 1.0 - 10 Hz,

Head cooling from -21 °C to 3 °C,

14" inch SMART TOUCH touch screen,

Menu in language: PL, ENG, GER, SWE, ESP, RUS,

Included:

Power cord,

Laser head,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for a period of 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

