



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **BECKMAN COULTER LITE Laboratory Centrifuge (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Beckman Coulter Microfuge-Lite  
is a 13,000 rpm table-top centrifuge equipped with an 18-place fixed-angle rotor,  
Specifications:  
Adjustable time from 1 to 15 minutes (continuous operation can be set),  
Maximum speed: 13,000 rpm,  
Includes 18 x 1.5 ml rotor,  
Includes:  
Power cord,  
BECKMAN F1802B Rotor,  
Dimensions: 9 x 11 x 7 cm,  
Weight: 4 KG,  
Has a current inspection and is ready for use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

