



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HETTICH EBA 20 laboratory centrifuge (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Hettich EBA 20
is a compact centrifuge suitable for small volume samples,
Technical specifications:
Rotor equipped with 8 cups of Ø 1.8 x 9 cm,
Maximum capacity: 8 x 15 mL,
Acceptable tube height: 13.5 cm,
Locking lid protects the centrifuge from accidental opening,
A stick is attached to the base of the centrifuge for emergency lid release in case of power failure,
Digital LED displays for speed, centrifugation time and centrifuge operation,
Foiled control panel,
Keypad on the control panel in English,
Speed control in the range: 500 - 6000 revolutions / min (step of 100 revolutions),
Centrifugation time range: 1 - 99 min (1 min step),
Possibility of short (< 1 min) centrifugation thanks to a button that activates the centrifuge for the time it is held down,
Noise level: ≤ 54 dB,
Included:
User manual in English, German, French, Italian (PDF),
Power cord,
Overall dimensions: 30 x 21 x 23 cm,
Weight: 4 kg,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

