



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**Magnetotherapy ASTAR MAGNER PLUS (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Polish production, (2013),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 500 VA,  
Magner Plus  
is a two-channel device used to perform magnetotherapy treatments. The channels are completely independent and can generate different fields at the same time,  
The user can use a bank of built-in programs or set all the parameters independently. You can also define and store your own treatment programs in the controller's memory,  
Application:  
Delayed bone fusion,  
Pseudarthritic joint,  
Osteoporosis,  
Joint degeneration,  
Rheumatoid arthritis,  
Ulcers and trophic changes of the lower extremities,  
Bacterial infections of the skin and soft tissues,  
Scarring,  
Post-stroke conditions,  
Migraine, vasomotor headache,  
Cranial and peripheral nerve dysfunction,  
Multiple sclerosis,  
Infectious disease of the cornea,  
Disappearance of the optic nerve,  
Hypertension,  
Coronary artery disease,  
Cardiac arrhythmia,  
Irritable colon,  
Long-term pancreatitis,  
Specifications:  
Maximum induction: 10mT/100Gs,  
Modulated mode parameters: pulse 1s/ interval 0.5-8s,  
Frequency range: 2-60Hz,  
Treatment timer: 1-30 min,  
MAGNETIC FIELD:  
Field types: sine, triangle, half-sine rectangle, half-triangle, half-rectangle,  
Mixed mode: different field types in one treatment session,  
Continuous and modulated emission,  
High magnetic field induction,  
Wide frequency range,  
MANUAL MODE:  
Full control over treatment parameters for advanced users,  
TREATMENT PROGRAMS:  
Simplifying the operation of the device,  
Disease units selectable by name,  
46 disease units,  
79 treatment programs,  
50 user-settable programs,  
ERGONOMY:  
Operating modes: program / manual,  
Character display - easy to read operation,  
Selection of disease units by name,  
Treatment timer,  
Two completely independent application channels,  
Four universal applicator slots,  
Automatic recognition of applicator type  
,  
CONSTRUCTION:  
Microprocessor control of camera operation,  
Liquid crystal character display,  
Durable and reliable keyboard,  
Autotest - ongoing control of the efficiency of the camera,  
Included:  
Power cord,  
PENTLA CSP35 spool applicator,  
Dimensions:  
Controller: 34 x 45 x 13.5 cm,  
Spool applicator: 35 cm diameter,  
Total weight: 20 KG,  
It has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

