



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HETTICH ROTOFIX 32A laboratory centrifuge (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Robust, and versatile centrifuge
HETTICH ROTOFIX 32A
for routine laboratory tasks in doctors' offices and small hospitals, clinical chemistry and cytology,
Rugged construction allows it to be used in industrial and research laboratories, for example, for testing chromium baths,
Features:
Metal housing and stainless steel rotating chamber,
Automatic rotor recognition,
Rotor unbalance protection,
Short spin button,
Ability to set time and speed and preview acceleration,
Specifications:
Mass capacity 4x 100ml / 32x 15ml,
Maximum speed: 6000 RPM,
Maximum acceleration: 4226 RCF,
Running time: 1-99 min, continuous operation, short spin button,
Minimum acceleration/braking time to maximum speed: 13/15,
Included:
Instruction manual in GER, ENG, FR, IT, PL, LT, LV, ET (PDF),
Power cord,
Tilt rotor
HETTICH 1624
,
4 pieces of adapters
HETTICH 1381
with cup inserts
HETTICH 1383
,
Dimensions: 27 x 43 x 26 cm,
Weight: 23 KG,
Cat 01,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

