



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



X-ray chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
X-ray stool

is designed for comfortable and stable positioning of the patient during X-ray imaging. The round upholstered seat and adjustable backrest and headrest provide ergonomics and comfort in any position. This X-ray chair is ideal for use in radiology labs,
Adjustable seat, backrest and headrest heights provide comfort so you can easily adjust the patient to a comfortable position and equipment requirements. The seat is mounted on a stable aluminum base with chrome actuator,

Features:

Round seat - comfort and easy cleaning,
Adjustable: seat, backrest and headrest height,
Chrome actuator - low or high variant,
Stable aluminum base,
Adjustable seat height in the range: 55 - 75 cm,
Free height adjustment of backrest and headrest with knobs,
Upholstery with medical certificate in green color (it is possible to replace the upholstery with a new one from our color pattern for an additional fee),

Dimensions:

Seat diameter: 40 cm,

Total height: 152 cm,

Weight: 8 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

