



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MPW MPW-350R Laboratory Centrifuge (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Polish production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 870 W,  
Centrifuge MPW-350R  
is used in medical, biochemical laboratories and food industry,  
Features of MPW-350R:  
Microprocessor control,  
Choice of 99 memory programs,  
Selection of 10 acceleration and deceleration curves,  
Automatic rotor identification,  
Error identification system,  
Registration of operating parameters - RS 232 (optional),  
Eco-friendly and efficient cooling system,  
Pre-cooling capability,  
Product compliance with international safety conditions,  
Technical specifications:  
RPM range: 300 - 15000 / rpm,  
Spinning time: 0 - 99 min,  
Temperature control range: -10 / + 40 °C,  
No rotor included - can be retrofitted at additional cost,  
Included:  
Power cord,  
Dimensions: 41 x 64 x 37 cm,  
Weight: 59 kg,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for a period of 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

