



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STRYKER 1115 PRIME SERIES transport cart (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Manufactured by U.S.A., (2016),
Transport cart
STRYKER PRIME SERIES 1115
is designed to reduce the physical burden on medical personnel, improve patient comfort and increase hospital efficiency,
New medical approved mattress in black,
Features:
Large wheels reduce the force required to move the wheelchair off the ground and maneuver,
Folding push handles,
Hydraulic system control available from three sides,
Side railing release function available at both ends of the wheelchair,
Hanger for IV fluids,
Technical specifications:
Height adjustment range in: 53 - 86 cm,
Back support in range: 0° - 90°,
Knee support: 0° - 40°,
Trendelenburg / Anti-Trendelenburg position: $\pm 17^\circ$,
Built-in central brake,
Dimensions:
Mattress: 62 x 200 cm,
Length of side rails: 36 x 147 cm,
Wheel diameter: 20 cm (8" inches),
Maximum load capacity: 318 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

